

Anti-TBCEL antibody (1-424aa) [AT1B10] (STJA0041728)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	219899
Gene Symbol	TBCEL
UniProt ID	TBCEL_HUMAN
Immunogen Region	1-424aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF staining of TBCEL in PC3 cells with STJA0041728 antibody Alexa Fluor 488 green and DAPI blue



Western blot of human TBCEL antibody STJA0041728 showing detection across multiple cell lysates



Flow cytometry of TBCEL in PC3 cells using Human TBCEL antibody with Alexa Fluor 488 detection