

Anti-PGAM1 antibody (1-254aa) [AT1G4] (STJA0041729)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2aL
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5223
Gene Symbol	PGAM1
UniProt ID	PGAM1_HUMAN
Immunogen Region	1-254aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of PGAM1 in HeLa cells stained with STJA0041729 antibody Alexa Fluor 488 secondary green and DAPI blue



Western blot of human PGAM1 protein expression across 293T Jurkat Raji A431 and HeLa cell lysates