

Anti-ACOT8 antibody (1-319aa) [AT4D10] (STJA0041730)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10005
Gene Symbol	ACOT8
UniProt ID	ACOT8_HUMAN
Immunogen Region	1-319aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ACOT8 in A549 cells with Alexa Fluor 488 green and DAPI blue



Western blot of Human ACOT8 protein using STJA0041730 antibody across multiple cell lysates



Flow cytometry of ACOT8 in CTLL-2 cells with goat anti-mouse IgG Alexa Fluor 488