

Anti-Fascin/FSCN1 antibody (1-493aa) [AT13D2] (STJA0041738)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6624
Gene Symbol	FSCN1
UniProt ID	FSCN1_HUMAN
Immunogen Region	1-493aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of FSCN1 in HeLa cells stained with STJA0041738 antibody. Alexa Fluor 488 secondary green and DAPI blue



Western blot of human FSCN1 Fascin protein in HeLa HepG2 and THP-1 cell lysates



Flow cytometry of FSCN1 in HeLa cells with human Fascin antibody STJA0041738