

Anti-FUBP1 antibody (279-448aa) [AT14F5] (STJA0041746)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8880
Gene Symbol	FUBP1
UniProt ID	FUBP1_HUMAN
Immunogen Region	279-448aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF staining of FUBP1 in HeLa cells with STJA0041746 antibody. Alexa Fluor 488 secondary green DAPI nuclear counterstain blue



Western blot of human FUBP1 protein in HeLa, HepG2 and Jurkat cell lysates