

## Anti-Glutathione S-transferase theta 1/GSTT1 antibody (1-240aa) [AT38D11] (STJA0041748)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">2952</a>
Gene Symbol	<a href="#">GSTT1</a>
UniProt ID	<a href="#">GSTT1_HUMAN</a>
Immunogen Region	1-240aa

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of GSTT1 in HepG2 cells stained with STJA0041748 antibody Alexa Fluor 488 secondary antibody and DAPI nuclear stain



Western blot of human GSTT1 antibody showing recombinant protein and cell lysates



Flow cytometry of GSTT1 in Hep3B cells stained with STJA0041748 antibody