

## Anti-PGK1 antibody (1-417aa) [AT2F4] (STJA0041756)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">5230</a>
Gene Symbol	<a href="#">PGK1</a>
UniProt ID	<a href="#">PGK1_HUMAN</a>
Immunogen Region	1-417aa

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of PGK1 in U87MG cells with DAPI and Alexa Fluor 488 staining



Western blot of human PGK1 antibody across HeLa LNCaP A549 HepG2 293T K562 and Jurkat cell lysates