

Anti-S100A4 antibody (1-101aa) [AT1C3] (STJA0041759)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6275
Gene Symbol	S100A4
UniProt ID	S10A4_HUMAN
Immunogen Region	1-101aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of S100A4 in A549 cells stained with STJA0041759 antibody Alexa Fluor 488 green and DAPI blue



Western blot of human S100A4 antibody using HeLa A549 NIH 3T3 and 3T3-L1 cell lysates



Western blot of human S100A4 antibody using HeLa A549 NIH 3T3 and 3T3-L1 cell lysates