

Anti-Glutathione S-transferase theta 2/GSTT2 antibody (1-244aa) [AT10H4] (STJA0041766)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2953
Gene Symbol	GSTT2
UniProt ID	GST2_HUMAN
Immunogen Region	1-244aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of GSTT2 in HeLa cells with DAPI and Alexa Fluor 488 staining



Western blot of human GSTT2 protein expression in multiple cell lines using STJA0041766 antibody



Flow cytometry of GSTT2 in HeLa cells with goat anti-mouse IgG Alexa Fluor 488 conjugate