

Anti-Nicotinamide N-Methyltransferase/NNMT antibody (1-264aa) [AT11G11] (STJA0041767)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4837
Gene Symbol	NNMT
UniProt ID	NNMT_HUMAN
Immunogen Region	1-264aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of NNMT in HeLa cells with DAPI nuclear staining and anti-NNMT antibody Alexa Fluor 488



Western blot of human NNMT antibody in HeLa and A549 cell lysates



Flow cytometry of NNMT in HeLa cells with goat anti-mouse IgG Alexa Fluor 488