

Anti-ALDOC antibody (1-364aa) [AT2E11] (STJA0041771)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	230
Gene Symbol	ALDOC
UniProt ID	ALDOC_HUMAN
Immunogen Region	1-364aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ALDOC in U87MG cells stained with STJA0041771 antibody



ICC IF of ALDOC in Jurkat cells with Alexa Fluor 488 and DAPI staining



Western blot of human ALDOC antibody showing PC3 U87MG WiDr Jurkat cell lysates and mouse brain tissue