

Anti-Peroxiredoxin 3/PRDX3 antibody (63-256aa) [AT1F8] (STJA0041781)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG3k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10935
Gene Symbol	PRDX3
UniProt ID	PRDX3_HUMAN
Immunogen Region	63-256aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of Peroxiredoxin 3 in HeLa cells with STJA0041781 antibody Alexa Fluor 488 secondary DAPI nuclear stain



Western blot of human PRDX3 protein using anti-PRDX3 antibody across multiple cell lines