

Anti-6-phosphogluconate Dehydrogenase/PGD antibody (1-483aa) [AT46B3] (STJA0041793)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5226
Gene Symbol	PGD
UniProt ID	6PGD_HUMAN
Immunogen Region	1-483aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of 6-phosphogluconate dehydrogenase in Hep3B cells with STJA0041793 antibody



Western blot of human PGD antibody STJA0041793 in HeLa, MCF7, 293T, A549, and Jurkat cell lysates