

Anti-Thioredoxin-2 antibody (60-166aa) [AT1B1] (STJA0041801)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	25828
Gene Symbol	TXN2
UniProt ID	THIOM_HUMAN
Immunogen Region	60-166aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of TXN2 in HeLa cells with DAPI and anti-human TXN2 antibody staining



Western blot of human Thioredoxin-2 protein expression across multiple cell lines



Western blot of human Thioredoxin-2 using STJA0041801 antibody with recombinant protein controls