

Anti-BAG2 antibody (1-211aa) [AT29E9] (STJA0041812)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9532
Gene Symbol	BAG2
UniProt ID	BAG2_HUMAN
Immunogen Region	1-211aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of BAG2 in HeLa cells stained with DAPI and Alexa Fluor 488



Western blot of Human BAG2 antibody STJA0041812 showing protein detection



Western blot of human BAG2 protein in Jurkat A549 K562 HepG2 and HeLa cell lysates