

Anti-BPGM antibody (1-259aa) [AT2E11] (STJA0041829)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	669
Gene Symbol	BPGM
UniProt ID	PMGE_HUMAN
Immunogen Region	1-259aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of BPGM in HeLa cells with Alexa Fluor 488 green and DAPI blue



Western blot of human BPGM antibody showing 293T and BPGM-transfected 293T cell lysates



Flow cytometry of BPGM in HeLa cells stained with STJA0041829 antibody