

Anti-ACADS antibody (25-412aa) [AT8B10] (STJA0041832)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	35
Gene Symbol	ACADS
UniProt ID	ACADS_HUMAN
Immunogen Region	25-412aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ACADS protein in Hep3B cells stained with STJA0041832 antibody Alexa Fluor 488 secondary green and DAPI nuclear stain blue



ICC IF of ACADS in A549 cells with STJA0041832 antibody Alexa Fluor 488 secondary DAPI nuclear stain



Western blot of human ACADS antibody in HepG2 and LnCap cell lysates