

Anti-MCM7 antibody (1-414aa) [AT118H7] (STJA0041854)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4176
Gene Symbol	MCM7
UniProt ID	MCM7_HUMAN
Immunogen Region	1-414aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of MCM7 in MCF7 cells stained with STJA0041854 antibody and Alexa Fluor 488



Western blot of human MCM7 antibody in HeLa Jurkat and K562 cell lysates



Western blot of human MCM7 antibody in HeLa Jurkat and K562 cell lysates