

Anti-TOMM20 antibody (25-145aa) [AT1B2] (STJA0041855)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9804
Gene Symbol	TOMM20
UniProt ID	TOM20_HUMAN
Immunogen Region	25-145aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of TOMM20 in HeLa cells with Alexa Fluor 488 and DAPI staining



Western blot of human TOMM20 protein expression across 9 cell lines using anti-human TOMM20 antibody



Flow cytometry of TOMM20 in HeLa cells using STJA0041855 antibody with Alexa Fluor 488 detection