

Anti-Sulfotransferase 2A1/SULT2A1 antibody (1-285aa) [AT13E10] (STJA0041864)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6822
Gene Symbol	SULT2A1
UniProt ID	ST2A1_HUMAN
Immunogen Region	1-285aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of SULT2A1 in Hep3B cells stained with STJA0041864 antibody



Western blot of human SULT2A1 protein in HepG2 cell lysate using anti-SULT2A1 antibody



Western blot of human SULT2A1 protein in 293T cell lysates using anti-SULT2A1 antibody