

Anti-Hexokinase 1/2/3 antibody (1-923aa) [AT4C12] (STJA0041868)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3101
Gene Symbol	HK3
UniProt ID	HXXK3_HUMAN
Immunogen Region	1-923aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of Human Hexokinase 1 2 3 antibody STJA0041868 across HeLa Jurkat K562 and mouse brain lysates



Western blot of Human Hexokinase 1 2 3 antibody STJA0041868 across HeLa Jurkat K562 and mouse brain lysates