

## Anti-Glutathione S-transferase pi 1/GSTP1 antibody (1-210aa) [AT12C10] (STJA0041871)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">2950</a>
Gene Symbol	<a href="#">GSTP1</a>
UniProt ID	<a href="#">GSTP1_HUMAN</a>
Immunogen Region	1-210aa

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of GSTP1 in A549 cells stained with STJA0041871 antibody Alexa Fluor 488 secondary and DAPI



Western blot of GSTP1 protein expression in K562 Jurkat HeLa 293T and A549 cell lysates



Flow cytometry of GSTP1 in K562 cells using STJA0041871 antibody with Alexa Fluor 488 detection