

Anti-XPA antibody (1-273aa) [AT71H3] (STJA0041875)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7507
Gene Symbol	XPA
UniProt ID	XPA_HUMAN
Immunogen Region	1-273aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of XPA in HeLa cells with STJA0041875 antibody Alexa Fluor 488 secondary DAPI nuclear stain



Western blot of human XPA antibody across SW480 K562 MCF7 LNCaP and Raji cell lysates