

Anti-Calcineurin A/PPP3CA antibody (1-511aa) [AT1E11] (STJA0041878)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5530
Gene Symbol	PPP3CA
UniProt ID	PP2BA_HUMAN
Immunogen Region	1-511aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of PPP3CA in HeLa cells stained with Human Calcineurin A antibody Alexa Fluor 488 green and DAPI blue



Western blot of human Calcineurin A PPP3CA antibody in HeLa 293T and LNCaP cell lysates