

Anti-IRF5 antibody (176-240aa) [10T1] (STJA0041891)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3663
Gene Symbol	IRF5
UniProt ID	IRF5_HUMAN
Immunogen Region	176-240aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF staining of IRF5 in THP-1 cells with STJA0041891 antibody Alexa Fluor 488 green and DAPI blue



ICC IF of IRF-5 in Raw264.7 cells stained with STJA0041891 antibody and Alexa Fluor 488



Western blot of human IRF5 antibody detecting IRF-5 protein in THP-1 cell lysate



Flow cytometry of IRF5 in THP-1 cells using STJA0041891 antibody with Alexa Fluor 488 detection