

Anti-BAK antibody (29-187aa) [AT38E2] (STJA0041896)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION



Gene ID	578
Gene Symbol	BAK1
UniProt ID	BAK_HUMAN
Immunogen Region	29-187aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of BAK in HeLa cells stained with STJA0041896 antibody Alexa Fluor 488 secondary green and DAPI nuclear stain blue



Flow cytometry of BAK protein in HeLa cells using STJA0041896 antibody with Alexa Fluor 488 detection