

## Anti-ISG15/UCRP antibody (1-157aa) [3E5] (STJA0041916)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/WB/ICC/IF    |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                          |
| Concentration       | 1 mg/mL                                                                                                                                             |
| Conjugation         | Unconjugated                                                                                                                                        |
| Purification        | By protein-A affinity chromatography                                                                                                                |
| Formulation         | Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol                                                                  |
| Isotype             | IgG1k                                                                                                                                               |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION


|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">9636</a>        |
| Gene Symbol      | <a href="#">ISG15</a>       |
| UniProt ID       | <a href="#">ISG15_HUMAN</a> |
| Immunogen Region | 1-157aa                     |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of ISG15 in HeLa cells stained with STJA0041916 antibody Alexa Fluor secondary antibody and DAPI nuclear stain



Western blot of human ISG15 UCRP protein in 488 MCF7 and PC3 cell lysates