

Anti-His-tag antibody [AT1D6] (STJA0041917)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	His Tag

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of His-tag antibody detecting OTUB2-His and GAD-His in transfected 293T cell lysates



Western blot of His-tag antibody detecting recombinant human TREM1 CD34 and ATG3 proteins