

Anti-GAPDH antibody [AT4E8] (STJA0041932)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Rabbit

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of GAPDH in HeLa cells stained with STJA0041932 antibody (Alexa Fluor 488 secondary) and DAPI nuclear stain



Western blot of GAPDH antibody across HepG2 MCF7 Jurkat K562 A431 LNCaP and Ramos cell lysates



Western blot of GAPDH antibody STJA0041932 in mouse brain spleen eye and muscle tissue lysates



Western blot of GAPDH antibody on 293T cell lysate with serial dilutions