

Anti-IRF3 antibody (108-166aa) [3F10] (STJA0041933)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3661
Gene Symbol	IRF3
UniProt ID	IRF3_HUMAN
Immunogen Region	108-166aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of IRF3 in HeLa cells stained with STJA0041933 antibody Alexa Fluor 488 green and DAPI blue



Flow cytometry of IRF3 in MCF7 cells using STJA0041933 antibody with Alexa Fluor 488 detection



Flow cytometry of IRF3 in Jurkat cells using STJA0041933 antibody with Alexa Fluor 488 detection



Flow cytometry of human IRF3 antibody STJA0041933 showing cell population distribution