

Anti-Hexokinase antibody (1-917aa) [4D7] (STJA0041942)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3098
Gene Symbol	HK1
UniProt ID	HXX1_HUMAN
Immunogen Region	1-917aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of Hexokinase in HeLa cells with Alexa Fluor 488 and DAPI staining



Western blot of recombinant human Hexokinase proteins 1-4 using anti-human Hexokinase antibody



Western blot of human Hexokinase antibody STJA0041942 across multiple cell lysates