

Anti-RalA/RalB antibody (1-203aa) [AT41E2] (STJA0041963)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5898
Gene Symbol	RALA
UniProt ID	RALA_HUMAN
Immunogen Region	1-203aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of recombinant human RalA and RalB proteins using anti-human RalA/RalB antibody



Western blot of human RalA/RalB antibody showing protein bands in HeLa, HepG2, NIH3T3, PC3 and U87MG cell lysates



Flow cytometry of RalA/RalB in HeLa cells using STJA0041963 antibody with Alexa Fluor 488 detection