

Anti-AGXT antibody (330-392aa) [AT2D10] (STJA0041977)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-A affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	189
Gene Symbol	AGXT
UniProt ID	AGT1_HUMAN
Immunogen Region	330-392aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC IF of AGXT in HeLa cells stained with STJA0041977 antibody Alexa Fluor 488 secondary antibody and DAPI nuclear stain



ICC IF of AGXT in Hep3B cells stained with STJA0041977 antibody Alexa Fluor 488 secondary green and DAPI nuclear stain blue



Western blot of recombinant human AGXT protein detected with anti-AGXT antibody



Flow cytometry of AGXT in HeLa cells using Human AGXT antibody STJA0041977