

Anti-Cellubrevin/VAMP-3 antibody (1-77aa) [AT4G9] (STJA0041988)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2bk
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9341
Gene Symbol	VAMP3
UniProt ID	VAMP3_HUMAN
Immunogen Region	1-77aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human cellubrevin VAMP-3 in mouse brain tissue lysate using STJA0041988 antibody



Western blot of human cellubrevin VAMP-3 in mouse brain tissue lysate using STJA0041988 antibody



Flow cytometry of Cellubrevin VAMP-3 in HeLa cells using STJA0041988 antibody with Alexa Fluor 488 detection