

Anti-visfatin/NAMPT antibody (1-491aa) [4D5] (STJA0041989)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG1k
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10135
Gene Symbol	NAMPT
UniProt ID	NAMPT_HUMAN
Immunogen Region	1-491aa

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC IF of visfatin in HeLa cells with STJA0041989 antibody. Alexa Fluor 488 secondary green and DAPI nuclear stain blue.



Western blot of human visfatin NAMPT in 293T cell lysate using anti-human visfatin antibody.



Western blot of human visfatin NAMPT protein in SW480 cell lysate using STJA0041989 antibody.