

Anti-Wnt-3a antibody (19-352aa) [3A6] (STJA0041990)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	By protein-G affinity chromatography
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol
Isotype	IgG2ak
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	89780
Gene Symbol	WNT3A
UniProt ID	WNT3A_HUMAN
Immunogen Region	19-352aa

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot of human Wnt-3a antibody detection in 293T cell lysates



Immunohistochemistry of human Wnt-3a antibody in breast cancer tissue with DAB detection