

Anti-B3GAT1 antibody [IHC539] (STJA0041993)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	1:100-1:200
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Gene ID	27087
Gene Symbol	B3GAT1
UniProt ID	B3GA1_HUMAN
Specificity	Tested and validated on tonsil tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---