

## Anti-Cytokeratin 7 antibody [IHC007] (STJA0042038)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	1:50-1:200
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgG2ak (79)
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">3855</a>
Gene Symbol	<a href="#">KRT7</a>
UniProt ID	<a href="#">K2C7_HUMAN</a>
Specificity	Tested and validated on salivary gland, lung adenocarcinoma tissue.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---