

## Anti-GLUT1 antibody [IHC404] (STJA0042043)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-p              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Clonality           | Monoclonal                                                  |
| Conjugation         | Unconjugated                                                |
| Dilution Range      | 1:100-1:400                                                 |
| Formulation         | Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide |
| Isotype             | IgG2a                                                       |
| Storage Instruction | Store at 2-8°C for up to 1 year. Do not freeze.             |

### TARGET INFORMATION

|             |                                                                              |
|-------------|------------------------------------------------------------------------------|
| Gene ID     | <a href="#">6513</a>                                                         |
| Gene Symbol | <a href="#">SLC2A1</a>                                                       |
| UniProt ID  | <a href="#">GTR1_HUMAN</a>                                                   |
| Specificity | Tested and validated on colorectal carcinoma, malignant mesothelioma tissue. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|