

Anti-SPNS2 antibody [IHC694] (STJA0042061)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	1:50-1:200
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Gene ID	124976
Gene Symbol	SPNS2
UniProt ID	SPNS2_HUMAN
Specificity	Tested and validated on intestine tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---