

## Anti-Calretinin antibody [IHC523] (STJA0042100)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-p              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Formulation         | Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide |
| Isotype             | IgG1  |
| Storage Instruction | Store at 2-8°C for up to 1 year. Do not freeze.             |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Gene ID     | <a href="#">794</a>                                    |
| Gene Symbol | <a href="#">CALB2</a>                                  |
| UniProt ID  | <a href="#">CALB2_HUMAN</a>                            |
| Specificity | Tested and validated on mesothelioma, appendix tissue. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|