

## Anti-CD34 antibody [IHC034] (STJA0042116)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-p              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Formulation         | Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide |
| Isotype             | IgG1  |
| Storage Instruction | Store at 2-8°C for up to 1 year. Do not freeze.             |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Gene ID     | <a href="#">947</a>                              |
| Gene Symbol | <a href="#">CD34</a>                             |
| UniProt ID  | <a href="#">CD34_HUMAN</a>                       |
| Specificity | Tested and validated on tonsil, placenta tissue. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|