

Anti-CD99 antibody [IHC126] (STJA0042131)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-p |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Formulation | Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide |
| Storage Instruction | Store at 2-8°C for up to 1 year. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|--|
| Gene ID | 4267 |
| Gene Symbol | CD99 |
| UniProt ID | CD99_HUMAN |
| Specificity | Tested and validated on ewing's sarcoma, pancreas, esophagus tissue. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|