

## Anti-CDK4 antibody [IHC077] (STJA0042132)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">1019</a>
Gene Symbol	<a href="#">CDK4</a>
UniProt ID	<a href="#">CDK4_HUMAN</a>
Specificity	Tested and validated on liposarcoma tissue.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---