

Anti-Cytokeratin 8/18 antibody [IHC559] (STJA0042156)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Rabbit/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	Rabbit IgG/Mouse IgG1
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Specificity	Tested and validated on breast ductal carcinoma, colon cancer tissue.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---