

## Anti-FOXP3 antibody [IHC773] (STJA0042172)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">50943</a>
Gene Symbol	<a href="#">FOXP3</a>
UniProt ID	<a href="#">FOXP3_HUMAN</a>
Specificity	Tested and validated on tonsil, lymph node tissue.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---