

Anti-Ki-67 antibody [IHC167] (STJA0042199)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgGk (Rab)
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Specificity	Tested and validated on tonsil tissue.
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---