

Anti-Nestin antibody [IHC105] (STJA0042224)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Gene ID	10763
Gene Symbol	NES
UniProt ID	NEST_HUMAN
Specificity	Tested and validated on tonsil, glioma tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---