

Anti-PD-L1 antibody [IHC441] (STJA0042241)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
UniProt ID	PD1L1_HUMAN
Specificity	Tested and validated on tonsil, lung adenocarcinoma tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---