

## Anti-PLAP antibody [IHC088] (STJA0042247)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-p              |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Clonality           | Monoclonal                                                  |
| Conjugation         | Unconjugated                                                |
| Formulation         | Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide |
| Storage Instruction | Store at 2-8°C for up to 1 year. Do not freeze.             |

### TARGET INFORMATION

|             |                                          |
|-------------|------------------------------------------|
| Gene ID     | <a href="#">9373</a>                     |
| Gene Symbol | <a href="#">PLAA</a>                     |
| UniProt ID  | <a href="#">PLAP_HUMAN</a>               |
| Specificity | Tested and validated on placenta tissue. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|